



DEPARTMENT OF HEALTH
Republic of South Africa

TEST PROCEDURES FOR FILM PROCESSING AND INTENSIFYING SCREENS

DEPARTMENT OF HEALTH

RADIATION CONTROL

<http://www.doh.gov.za/department/radiation/01.html>

Contents

1	DARKROOM SAFELIGHT TEST (FOG TEST)	3
2	HOW TO CALCULATE SENSITOMETRIC CONTROL PARAMETERS	5
3	PROCEDURE FOR ESTABLISHMENT OF A SENSITOMETRIC PROCESSOR CONTROL PROGRAMME	7
4	PROCEDURE FOR DAILY SENSITOMETRIC CONTROL OF SPEED, CONTRAST AND BASIC FOG FOR AUTOMATIC PROCESSOR	10
5	CHANGING OF CONTROL EMULSIONS (FILM CROSSOVER PROCEDURES)	14
6	TEST FOR RESIDUAL THIOSULPHATE	15
7	VISUAL INSPECTION, IDENTIFICATION AND CLEANING OF CASSETTES	16
8	ASSESSMENT OF THE LIGHT TIGHTNESS OF CASSETTES	17
9	VISUAL INSPECTION, IDENTIFICATION AND CLEANING OF INTENSIFYING SCREENS	18
10	SATISFACTORY CONTACT BETWEEN X-RAY FILM AND INTENSIFYING SCREENS	19
11	ASSESSMENT OF THE RELATIVE SPEED OF INTENSIFYING SCREENS FOR REPLACEMENT	20
12	ASSESSMENT OF THE RELATIVE SPEED OF INTENSIFYING SCREENS IN CURRENT USE	21
13	ASSESSMENT OF THE RELATIVE SPEED OF INTENSIFYING SCREENS WHEN NEW SCREENS ARE INTRODUCED	22
14	FORMS	23
14.1	PROCESSING MAINTENANCE STANDARDS	23
14.2	MAINTENANCE STANDARDS DENSITY READINGS	24
14.3	CROSSOVER WORKSHEET	25
14.4	CROSSOVER WORKSHEET - EXAMPLE	26
14.5	DAILY PROCESSOR CONTROL CHART:	27
14.6	X-RAY PROCESSING CONTROL CHART ACTION FORM	28
14.7	INDIVIDUAL CASSETTE AND SCREEN RECORD	29
14.8	RECORD OF QC TESTS: CASSETTES AND SCREENS	30
14.9	RECORD OF QC TESTS: CASSETTES AND SCREENS	31
14.10	REPEAT AND REJECT ANALYSIS	32
14.11	SCREEN SPEED EVALUATION TEST RESULTS	33

1 DARKROOM SAFELIGHT TEST (FOG TEST)

PURPOSE

To ensure that the conditions in the darkroom complies with the minimum conditions as prescribed

EQUIPMENT REQUIRED

Loaded cassette (18 x 24 cm)

Densitometer

Cardboard

Stopwatch / Timer

Tape measure

X-Ray unit

PROCEDURE

1. Expose the loaded cassette to radiation to produce an optical density of $\pm 0.8 - 1.0$.
 - ± 40 kVp and 2 mAs at 100 cm FFD for a 400-speed system to sensitize the film emulsion.
2. In the darkroom, switch off all safelights.
3. Unload the film from the cassette and place it on the work bench or area to be assessed.
4. Cover one half of the film with the piece of cardboard.
5. Switch on the safelight(s), activate the stopwatch and wait for a period of two minutes.
6. Process the film.

IMAGE

The film will display two areas of uniform densities, divided by a straight line (line represented by the edge of the cardboard that had covered the one half of the film).

INTERPRETATION

1. Zero and calibrate the densitometer.
2. Measure the densities on each half of the film (± 1 cm from the dividing line).
3. The measured density difference may not exceed a maximum of 0.05.

CORRECTIVE ACTION

1. Verify the distance from the safelight(s) to the work bench or feeding tray(s).
 - Distance must be at least 120 cm.

2. Check the wattage of light bulbs in the safelights(15 – 25 Watts).
3. Verify that the correct safelight filters are used for the film types in use.
 - U/V-Blue Film (Monochromatic): *Amber Filter*(mixture of red, orange and some yellow).
 - U/V-Blue-Green (Orthochromatic): *Pure Red Filter*.
 - Multifomat- or Laser camera Film (Red or amber sensitivity): *Green Filter*.
 - Panchromatic Film (Cine-Film): *Total darkness*.
4. Inspect the safelight filters for cracks or pinholes.

A process of elimination must be used to correct each of the above parameters and the test must be repeated until the requirements are met.

2 HOW TO CALCULATE SENSITOMETRIC CONTROL PARAMETERS

PROCEDURE

1. Expose a film on both sides using a sensitometer. (Only emulsion side for single emulsion film – two exposures).
2. Process the sensitometric control strip.
3. Measure the densities using a densitometer.

CONTRAST INDEX

1. Select two steps within the useful density range (0.25 above base + fog & 2.0 above base + fog).
2. The steps should not be adjacent but should show an adequate density difference, e.g. 4 - 7 step separation.
3. Measure and record the densities of the two steps and calculate the density difference between these steps.
4. The same two steps should be monitored daily.
5. The contrast index is plotted daily on a graph.

Acceptable variation: Density values of 0.15 above or below the speed index.

SPEED INDEX

1. Measure the densities on the strip with a densitometer.
2. Select a step with a density of 1 above base + fog (towards the middle of the straight line portion of the characteristic curve) under normal circumstances.
3. Measure and record the density of this step.
4. This selected step will be the step that is measured daily to monitor the speed index.
5. Daily density readings of this step are plotted on a graph.

Acceptable variation: Density values of 0.15 above or below the speed index.

BASE AND FOG LEVEL

1. The density is measured in the part of the film, which was shielded from intentional exposure when the test strip was made (e.g. centre of the film or step 1).
2. Plot this value (fog index) on a graph.

Any change taken place will be towards an increase in fog level. A fog level above 0.3 is unacceptable.

If the fog level enters the useful density range there will be deterioration in contrast.

Please note that once established, these index values should not change unless a deliberate change or changes have been made to the processing standards, for example:

- the introduction of a new film type;
- the use of chemistry from a different supplier;
- wilful adjustment of the developer temperature;
- wilful adjustment of the developer and/or fixer replenishment rates, and
- change in the processing cycle time.

3 PROCEDURE FOR ESTABLISHMENT OF A SENSITOMETRIC PROCESSOR CONTROL PROGRAMME

PURPOSE

To determine the standard operating levels for all processors in the department – i.e. STANDARD PROCESSING CONDITIONS. These levels should be the same for all processors processing the same film type.

EQUIPMENT REQUIRED

Sensitometer
Densitometer
Digital thermometer
Stop watch
Control emulsion (film type most frequently used)
Processor and fresh chemicals
Control chart

PREPARATION

Darkroom

Darkroom must be checked for correct safelight conditions and white light leaks.

Processor

Perform all the actions of a processor service according to the manual of your processor and the recommendations of the supplier of your chemicals, which should include the following:

1. Clean the processor rollers and replenishment tanks thoroughly. Flush out all lines, including the tubing from the replenishment.
2. Replace the developer re-circulation filter.
3. Use a checklist to inspect all components of the processor. Inspect alignment and wear points.
4. Correct any mechanical problems. Make note of and order new parts when indicated. Record this information.
5. Lubricate all necessary points.
6. Carefully mix and install fresh chemicals in the following order:
 - a. Fixer and developer replenishment tanks. While operating the replenisher pumps temporarily to assure that all fresh water is flushed out of the replenisher lines and to assure that the replenisher pumps are functioning properly.
 - b. Flush processor fixer tank with fresh water.

- c. Fill fixer tank and replace rack.
 - d. Flush developer tank with fresh water.
 - e. Fill developer tank and replace rack. Make note of chemical brand, type, order code or other identification and serial number. **Also note amount of starter solution added to the developer solution in the tank.**
7. Replace remaining racks.
 8. Perform start-up procedure.
 9. Allow the processor to operate for 30 minutes.
 10. Check and record:
 - a. Developer temperature
 - b. Fixer temperature
 - c. Wash water temperature. The chemistry temperatures should be within 0.5 °C of those recommended
 - d. Replenishment rates
 - e. Film transport time
 11. Process films and wait until approximately 50 (35 x 43 cm) films have been processed or the developer is seasoned.
 - a. To season the developer, use "starter" solution obtainable from the chemical supplier and add the required quantity to the developer solution in the developer processing tank (amount in millilitre per developer tank volume of your processor) to season the developer.

PROCEDURE

1. Expose at least six sensitometric films with the sensitometer under the recommended safelight conditions. Double emulsion film should be exposed on both sides of the film. Single emulsion film should be exposed only on the emulsion side.
2. Wait at least 30 minutes but not longer than 4 hours before processing of the films.
3. Process the films - crosswise, same position on feed tray (this position should be standardised).
4. Zero and calibrate the densitometer.
5. Read the densities on the six strips. Be sure to read the densities in the centre of each strip, not near the edges. Check the zero and calibration of the densitometer after reading each strip.
6. Determine the average of the densities of the six strips.
7. Select and mark the steps producing the densities nearest to 0.25; 1.0 and 2.0 above the base-plus-fog level.

8. Determine the appropriate values and record these values and indicate the control limits on the control chart.

Three values should be recorded:

- Contrast Index (± 0.15)
 - Speed Index (± 0.15)
 - Base and Fog (± 0.03) and less than 0.3
9. Record all operating levels on the processor maintenance standard form and the control chart (graph).
 10. Adjust each processor in the department that uses the same type of chemistry and film to produce the same density for the same steps on the sensitometric control strip that have received the same exposure.

CORRECTIVE ACTION

If you have difficulties with any of the above procedures, contact the technical service representative from the firm that manufactures your film and chemistry.

4 PROCEDURE FOR DAILY SENSITOMETRIC CONTROL OF SPEED, CONTRAST AND BASIC FOG FOR AUTOMATIC PROCESSOR

PURPOSE

Sensitometric control would indicate variations in the processing conditions before they are noticed on the processed radiographs and the cause for the variation can thus be rectified before marked changes are visible on the radiographs or before breakdown of the processing cycle occurs.

REQUIREMENTS

Sensitometer
 Densitometer
 Digital thermometer
 Control emulsion

CONDITIONS AND DATA TO BE RECORDED

1. The standardised processing conditions previously established must be adhered to.
2. The following data must be recorded on:
 - a. The sensitometric film
 - Date and time film was processed
 - Colour of sensitometric exposure
 - Processor being evaluated
 - Measured Speed Index, Contrast Index, Base and Fog level
 - b. The control chart
 - Speed Index, Contrast Index, Base and Fog.
 - This information is charted in terms of **deviation from** established normal values. For example, if the normal *Speed Index* density = 1.15 and the sensitometric reading = 1.20 it would be charted as +0.05 from the normal value.
 - Processor Maintenance Standards
 - Developer temperature
 - Fixer temperature
 - Wash water temperature
 - Developer Replenishment rate
 - Fixer Replenishment rate

- Dryer temperature

CONTROL INDEXES

• SPEED INDEX

- Speed index refers to that density step chosen in the standardization procedure which has a density of ± 1 above base + fog.
- Acceptable variation is: 15 - 20 % i.e. 0.15 above or below the standardised speed index.

• CONTRAST INDEX

- Contrast index refers to the density differences between two steps, in the useful density range (0.25 - 2.0 above base + fog), chosen in the standardization procedure. The steps should not be adjacent but should show an adequate density difference (4 - 7 steps).
- Acceptable variation is: 15 - 20 % i.e. 0.15 above or below the standardized contrast index.

• BASE AND FOG

- Base and fog refers to the density measured in part of the film, which was shielded from intentional exposure when the test strip was made.
- An expected change will be towards an increase in the fog level.
- A fog level above 0.25 is suspect and a fog level of 0.30 is unacceptable.

METHOD

1. Follow the manufacturer's start-up procedure everyday.
2. Run several clean-up sheets and check them for roller marks and scratches (Use exposed unprocessed film).
3. Allow sufficient time for the processor temperature to stabilize - $\pm \frac{1}{2}$ hour.
4. Check the following:
 - a. Solution temperatures
 - b. Replenishment rates
 - c. Water flow rates
 - d. Dryer temperature
5. Expose the control film(s) with the sensitometer under safelight conditions. Dual emulsion film should be exposed on both emulsions.
6. Process the film (Feed film at the standardized position on the feed tray e.g. crosswise / lengthwise and left / centre / middle).

IMAGE

The processed film will show two series of stepped exposures.

INTERPRETATION

1. Zero and calibrate the densitometer.
2. Read the density of the base-fog level.
3. Read the density levels of the two steps representing the contrast index on each strip.
4. Read the density level of the step representing the speed index on each strip.
5. Average the values for the two strips on the film.
6. Plot values on the control chart.
7. Analyse the control charts carefully.
 - a. Are all three points within the control limits?
 - b. Are there any apparent trends?
8. If any single point (or points) falls outside of the control limits, run two more strips and verify that the first readings were correct. If points are still outside the limits, corrective action must be taken before any patient films are processed.
9. Record relevant data on the control charts.

CORRECTIVE ACTION

1. Make the adjustments that you believe will bring the processor back into control.
 - a. If control limits are exceeded check the following:
 - i. Developer temperature
 - ii. Replenishment rates
 - iii. Water flow rates
 - iv. Water temperature (if appropriate)
 - v. Re-circulation
 - vi. Filters
 - vii. Batch mix dates
 - viii. Recent maintenance
 - ix. Film fog
 - x. Transport time
 - xi. Control emulsion
 - b. If trends are noted in the control chart, check the following:

- i. Developer temperature
 - ii. Replenishment rates
 - iii. Change in mix, types and number of films processed
 - iv. Proper mixing of replenisher
 - v. Control emulsion age or fog
 - vi. Leaks from overflow from the fixer tank getting into the developer tank
 - vii. Gremlins
2. Make only one adjustment to the processor at a time.
 3. After each change run another sensitometric strip.
 4. Record the types of changes made and the resultant change in the three control parameters.
 5. Special considerations (Miscellaneous)
 - A processor quality control programme can only be effective if the processors are maintained in a clean, functional condition as intended by the manufacturer.
 - A sensitometric control strip should be processed after any maintenance or cleaning operation.
 - All maintenance and cleaning should be recorded on a log that is maintained near the processor, since this may provide vital clues to otherwise unexplainable shifts in photographic processor quality.

5 CHANGING OF CONTROL EMULSIONS (FILM CROSSOVER PROCEDURES)

PURPOSE

When one nears the end of the box of film that has been used for processor monitoring, crossover testing should begin with the new box of film.

CHANGE OVER FROM ONE BOX TO THE NEXT (SAME BATCH/EMULSION NUMBER)

With the processor within the predetermined limits, one piece of film from the old box and one piece of film from the new box should be exposed sensitometrically and processed with the same orientation. Gross fog, speed and contrast should be read on both films, and the correction factor for the new film should be calculated to adjust the limits.

CHANGE OVER FROM ONE BOX TO THE NEXT (DIFFERENT BATCH/EMULSION NUMBER)

When changing from one batch of film to the other, six control strips of the old and six strips of the new emulsion should be processed.

The average values for the three basic parameters to be controlled should be determined and the standardised values be changed accordingly on the control chart.

6 TEST FOR RESIDUAL THIOSULPHATE

PURPOSE

To assess the efficiency of the washing process thus determining how effective residual soluble fixing salts are removed from the film. Salts remaining on the film will deteriorate the image in time.

REQUIREMENTS

Chemical Reagent Test-pack with comparator card

Processed radiograph

METHOD

1. Select a clear area on the film not carrying a developed density (clear film base).
2. Spot the reagent on the selected area by means of a dropper (one or two drops).
3. Leave for two minutes.
4. Blot the reagent off - a brownish stain will be left on the tested area.
5. **Immediately** compare the colour of the stain with the set of tints on the comparator chart in good, ordinary room lighting.

NB The stain will darken when exposed to light, which is bright, especially if it is direct sunlight. A long delay in assessing the depth of the colour will alter the result of the test.

INTERPRETATION

1. The card bears a set of tints, which begins with a pale straw colour (corresponding to a low level of residual thiosulphate) and ends with a light tan (corresponding to a higher level of residual thiosulphate).
2. Information on the comparator chart relates the tones shown on it to amounts of thiosulphate and the levels of acceptability for keeping radiographs for normal use and for archival storage.
3. A thiosulphate level of < 3 micro-g/ sq. cm. is recommended for archival permanence.
4. If the result of the test shows the washing to be inefficient, the fault may lie with one of the following:
 - a. water flow rate inadequate
 - b. water circulation pump malfunctioning
 - c. wash water too cold.
5. The level must be lower when the keeping-time is longer.
6. Radiographs kept as medical records in hospital fit into the group classified for ordinary use.

7 VISUAL INSPECTION, IDENTIFICATION AND CLEANING OF CASSETTES

PURPOSE OF TEST

To ensure the cassettes are in good order, are clearly marked with the correct type & speed of intensifying screens, and are identified by number to enable any cassette causing a film fault to be easily traced.

EQUIPMENT REQUIRED

Cassettes to be checked.

Intensifying screen identification labels (provided by manufacturer).

Indelible marker.

METHOD

1. Clean external surfaces of cassettes with damp swabs and soap or spirit-based swabs, taking particular care to avoid spirit contact with intensifying screens.
2. Inspect cassette for damage, particularly the edges, hinge and catches.
3. Check that the type of intensifying screen is clearly marked.
4. Check that the cassette is marked with a number, which corresponds to the number on the intensifying screen.

ASSESSMENT AND EVALUATION

If the cassette is damaged in a way that may allow light to enter the cassette it should be tested for light tightness. Cassettes with damage to the hinge or catches should be sent for repair or discarded. Intensifying screen identification labels should be replaced if not clear. All cassettes should be marked with indelible marker with a number to correspond to the number on the intensifying screen.

8 ASSESSMENT OF THE LIGHT TIGHTNESS OF CASSETTES

PURPOSE OF TEST

To test for light leakage into cassettes causing film fogging.

EQUIPMENT REQUIRED

Cassette to be tested.

Light source, e.g. x-ray film illuminator.

METHOD

1. Load the cassette with a film
2. Place the cassette with suspect area of cassette uppermost, next to light source.
3. Leave for 15 - 30 minutes.
4. Process the film.

ASSESSMENT AND EVALUATION

Any light leakage will cause an area of density on the film. If this light fogging is likely to interfere with an image when the cassette is in use, the cassette must be sent for repair or discarded.

9 VISUAL INSPECTION, IDENTIFICATION AND CLEANING OF INTENSIFYING SCREENS

PURPOSE OF TEST

To ensure that the intensifying screens are undamaged, are clearly identified, so that any artefact found on a radiograph can be traced to the offending cassette, and marks are removed.

EQUIPMENT REQUIRED

Cassette to be checked.

Method of marking screens, e.g. Letraset, biro.

Cleaning fluid recommend by intensifying screen manufacturer.

METHOD

1. Open cassette and inspect intensifying screen in bring light, looking for marks on the screen surface and any damage to the screen super coat. If the surface of the intensifying screen appears to be damaged a test radiograph must be taken to assess the extent of any artefact.
2. Check that the intensifying screens are marked with a number which corresponds to the number on the outside of the cassette.
3. Clean the intensifying screens using fluid recommended by manufacturer and following the instructions carefully.
4. Leave cassette open to dry for 30 minutes

ASSESSMENT AND EVALUATION

All the intensifying screens in use should be free from any artefacts and clearly marked for identification purposes. Discard any intensifying screens which cause artefacts. Any intensifying screens that are not identified should be marked using either "Letraset" or a biro. Note that this mark will be permanent, therefore great care must be taken in the marking process. This number should be positioned so as not to impede on the image but be visible within an exposed area. A suggested place is close to the identification window in the cassette.

10 SATISFACTORY CONTACT BETWEEN X-RAY FILM AND INTENSIFYING SCREENS

PURPOSE OF TEST

To assess any unsharpness in the image as a result of poor contact between the x-ray film and intensifying screens.

EQUIPMENT REQUIRED

Cassette to be tests.

Test object - wire mesh, size 36x44 cm, i.e. large enough to cover all cassettes. Wire diameter about 0.5 mm and mesh spacing about 6 times wire diameter x-ray generator and tube.

METHOD

1. Load the cassette to be tested with an x-ray film.
2. Place the test object on top of the cassette.
3. Position the x-ray tube so that the x-ray beam covers the cassette, with the FFD at least 1 m. A distance of 1.5 m is preferable as this minimizes geometrical unsharpness.
4. Expose the film. The exposure factors should be such that the kilovolt age does not penetrate the wire mesh (approximately 50 kVp) and the density of the image should be between 1.0 and 2.0. The correct density depends upon the closeness (gauge) of the mesh (see assessment).
5. Process the radiograph.

ASSESSMENT AND EVALUATION

Place the film on an x-ray film illuminator and view at normal viewing distance. The outline of the mesh should be clearly seen with a sharp contrast between the densities of the wire and the unfiltered beam. When a close mesh is used the background density needs to be darker (approximately 2.0) and any area of unsharpness will be clearly seen as an area of increased density. Many of the test objects available have amore open mesh and when these are used the density should be about 1.0. The area of unsharpness is then more easily seen as poorly defined edges to the image of the mesh. Any cassettes in which poor film-intensifying screen contact is demonstrated in an area likely to affect the image should be replaced.

11 ASSESSMENT OF THE RELATIVE SPEED OF INTENSIFYING SCREENS FOR REPLACEMENT

PURPOSE OF TEST

To assess any significant variation in intensifying screen speeds when new screens are being introduced into the department.

EQUIPMENT REQUIRED

Cassettes with intensifying screens in present use
Step wedge or phantom of uniform thickness
Cassette with intensifying screens to be evaluated
Densitometer (optional)

METHOD

1. Load both cassettes with film from the same box.
2. Place cassettes side by side on the tabletop.
3. Place the step wedge or phantom so that it covers part of each cassette.
4. Centralize the x-ray beam over the cassettes and collimate to cover the step wedge/phantom.
5. Set exposure factors.
6. Make an exposure.
7. Process both films in the same processor simultaneously if possible.

ASSESSMENT AND EVALUATION

The radiographs should be viewed side by side on an illuminator (viewing box). The image densities should not vary significantly when viewed with the naked eye. Densitometer readings should be within 10% of each other

Note that the attenuation of cassettes can vary from type to type. To obtain consistent results all the cassettes in the department containing intensifying screens of the same speed should be of the same type. Any variation in speed may be caused by a change of response in the intensifying screens in use. Repeat the test against a different cassette to assess the consistency of speed loss. Discuss the results with the manufacturer of the intensifying screens.

12 ASSESSMENT OF THE RELATIVE SPEED OF INTENSIFYING SCREENS IN CURRENT USE

PURPOSE OF TEST

To identify a difference in speed in intensifying screens of the same type in current use in the department.

EQUIPMENT REQUIRED


Cassettes with intensifying screens to be tested

Phantom of uniform thickness (at least 30 x 30 cm in area)

Densitometer (optional)

METHOD

1. Load all the cassettes with film from the same box.
2. Place four cassettes corner to corner on the tabletop.
3. Place the phantom in the centre of the four cassettes so that it covers one corner of each cassette.



4. Centralize the x-ray beam over the cassettes and collimate to cover the phantom.
5. Set exposure factors.
6. Make an exposure.
7. Repeat for all cassettes to be tested, keeping one of the cassettes unchanged throughout.
8. Process the films as soon as possible using the same processor.

ASSESSMENT AND EVALUATION

Providing the cassettes are all of the same type, the density produced on the radiograph should not vary significantly. Densitometer readings should be within 10% of each other. Any significant variation should be discussed with the manufacturer of the intensifying screens.

13 ASSESSMENT OF THE RELATIVE SPEED OF INTENSIFYING SCREENS WHEN NEW SCREENS ARE INTRODUCED

PURPOSE OF TEST

To assess any significant variation in intensifying screen speeds when new screens are being introduced into the department.

EQUIPMENT REQUIRED

Cassettes with intensifying screens in present use

Step wedge or phantom of uniform thickness

Cassette with intensifying screens to be evaluated

Densitometer (optional)

METHOD

1. Load both cassettes with film from the same box.
2. Place cassettes side by side on the tabletop.
3. Place the step wedge or phantom so that it covers part of each cassette.
4. Centralize the x-ray beam over the cassettes and collimate to cover the step wedge/phantom.
5. Set exposure factors.
6. Make an exposure.
7. Process both films in the same processor simultaneously if possible.

ASSESSMENT AND EVALUATION

The radiographs should be viewed side by side on an illuminator (viewing box). The image densities should not vary significantly when viewed with the naked eye. Densitometer readings should be within 10% of each other.

Note that the attenuation of cassettes can vary from type to type. To obtain consistent results all the cassettes in the department containing intensifying screens of the same speed should be of the same type. Any variation in speed may be caused by a change of response in the intensifying screens in use. Repeat the test against a different cassette to assess the consistency of speed loss. Discuss the results with the manufacturer of the intensifying screens.

14 FORMS

14.1 PROCESSING MAINTENANCE STANDARDS

PROCESSOR / DARKROOM:					PERIOD:				
CHEMICALS									
	Type			Supplier					
Developer									
Starter									
Fixer									
PROCESSOR									
	Vol. (L/ml)	Cycle Time (sec.)	Temp. (°C)	Repl. rate (ml/area or vol.)	Total Cycle Time (sec.)	Verification Date	Verification Date	Verification Date	Verification Date
Developer									
Starter									
Fixer									
Water									
Dryer									
SENSITOMETRIC CONTROL STANDARDS									
QA CONTROL FILMS									
Cross-over Date									
Type									
Batch no.									
Expiry date									
SENSITOMETRIC VALUES									
Parameters	Base +fog	Speed Index	Contrast Index	Base +fog	Speed Index	Contrast Index	Base +fog	Speed Index	Contrast Index
Date									
Step no.									
Density Value									

14.2 MAINTENANCE STANDARDS DENSITY READINGS

Hospital/institution:					Date:					
Processor:			Cycle time(s):			Temp.:				
Chemicals:										
Sensitometric Conditions:					e.g. Red Safelight; Green; X-Rite Dual Colour Model 334					
Film Feeding Orientation:					e.g. Emulsion up; lengthwise (notches first); Left side of feed tray					
Film Type(s):			Batch no.:			Expiry date:				
Sensitometer s/n:					Densitometer s/n:					

<table border="1" style="width:100%; border-collapse: collapse;"> <tr> <td style="width:10%;">Film no.</td> <td colspan="7">Density Values:</td> <td style="width:10%;"></td> <td style="width:10%;">Film no.</td> <td colspan="7">Density Values:</td> </tr> <tr> <td>1</td> <td>B + F =</td> <td></td> <td>Step no.</td> <td>Left</td> <td>Right</td> <td>Top</td> <td>Bottom</td> <td>Average</td> <td>4</td> <td>B + F =</td> <td></td> <td>Step no.</td> <td>Left</td> <td>Right</td> <td>Top</td> <td>Bottom</td> <td>Average</td> </tr> <tr> <td></td> <td colspan="7">Nearest to 1 + (B+F)</td> <td></td> <td></td> <td colspan="7">Nearest to 1 + (B+F)</td> </tr> <tr> <td></td> <td colspan="7">Nearest to 0.5</td> <td></td> <td></td> <td colspan="7">Nearest to 0.5</td> </tr> <tr> <td></td> <td colspan="7">Nearest to 2.0 + B+F</td> <td></td> <td></td> <td colspan="7">Nearest to 2.0 + B+F</td> </tr> <tr> <td></td> <td colspan="7">Contrast Index</td> <td></td> <td></td> <td colspan="7">Contrast Index</td> </tr> </table>	Film no.	Density Values:								Film no.	Density Values:							1	B + F =		Step no.	Left	Right	Top	Bottom	Average	4	B + F =		Step no.	Left	Right	Top	Bottom	Average		Nearest to 1 + (B+F)									Nearest to 1 + (B+F)								Nearest to 0.5									Nearest to 0.5								Nearest to 2.0 + B+F									Nearest to 2.0 + B+F								Contrast Index									Contrast Index							<table border="1" style="width:100%; border-collapse: collapse;"> <tr> <td style="width:10%;">Film no.</td> <td colspan="7">Density Values:</td> <td style="width:10%;"></td> <td style="width:10%;">Film no.</td> <td colspan="7">Density Values:</td> </tr> <tr> <td>2</td> <td>B + F =</td> <td></td> <td>Step no.</td> <td>Left</td> <td>Right</td> <td>Top</td> <td>Bottom</td> <td>Average</td> <td>5</td> <td>B + F =</td> <td></td> <td>Step no.</td> <td>Left</td> <td>Right</td> <td>Top</td> <td>Bottom</td> <td>Average</td> </tr> <tr> <td></td> <td colspan="7">Nearest to 1 + (B+F)</td> <td></td> <td></td> <td colspan="7">Nearest to 1 + (B+F)</td> </tr> <tr> <td></td> <td colspan="7">Nearest to 0.5</td> <td></td> <td></td> <td colspan="7">Nearest to 0.5</td> </tr> <tr> <td></td> <td colspan="7">Nearest to 2.0 + B+F</td> <td></td> <td></td> <td colspan="7">Nearest to 2.0 + B+F</td> </tr> <tr> <td></td> <td colspan="7">Contrast Index</td> <td></td> <td></td> <td colspan="7">Contrast Index</td> </tr> </table>	Film no.	Density Values:								Film no.	Density Values:							2	B + F =		Step no.	Left	Right	Top	Bottom	Average	5	B + F =		Step no.	Left	Right	Top	Bottom	Average		Nearest to 1 + (B+F)									Nearest to 1 + (B+F)								Nearest to 0.5									Nearest to 0.5								Nearest to 2.0 + B+F									Nearest to 2.0 + B+F								Contrast Index									Contrast Index						
Film no.	Density Values:								Film no.	Density Values:																																																																																																																																																																																																					
1	B + F =		Step no.	Left	Right	Top	Bottom	Average	4	B + F =		Step no.	Left	Right	Top	Bottom	Average																																																																																																																																																																																														
	Nearest to 1 + (B+F)									Nearest to 1 + (B+F)																																																																																																																																																																																																					
	Nearest to 0.5									Nearest to 0.5																																																																																																																																																																																																					
	Nearest to 2.0 + B+F									Nearest to 2.0 + B+F																																																																																																																																																																																																					
	Contrast Index									Contrast Index																																																																																																																																																																																																					
Film no.	Density Values:								Film no.	Density Values:																																																																																																																																																																																																					
2	B + F =		Step no.	Left	Right	Top	Bottom	Average	5	B + F =		Step no.	Left	Right	Top	Bottom	Average																																																																																																																																																																																														
	Nearest to 1 + (B+F)									Nearest to 1 + (B+F)																																																																																																																																																																																																					
	Nearest to 0.5									Nearest to 0.5																																																																																																																																																																																																					
	Nearest to 2.0 + B+F									Nearest to 2.0 + B+F																																																																																																																																																																																																					
	Contrast Index									Contrast Index																																																																																																																																																																																																					
<table border="1" style="width:100%; border-collapse: collapse;"> <tr> <td style="width:10%;">Film no.</td> <td colspan="7">Density Values:</td> <td style="width:10%;"></td> <td style="width:10%;">Film no.</td> <td colspan="7">Density Values:</td> </tr> <tr> <td>3</td> <td>B + F =</td> <td></td> <td>Step no.</td> <td>Left</td> <td>Right</td> <td>Top</td> <td>Bottom</td> <td>Average</td> <td>6</td> <td>B + F =</td> <td></td> <td>Step no.</td> <td>Left</td> <td>Right</td> <td>Top</td> <td>Bottom</td> <td>Average</td> </tr> <tr> <td></td> <td colspan="7">Nearest to 1 + (B+F)</td> <td></td> <td></td> <td colspan="7">Nearest to 1 + (B+F)</td> </tr> <tr> <td></td> <td colspan="7">Nearest to 0.5</td> <td></td> <td></td> <td colspan="7">Nearest to 0.5</td> </tr> <tr> <td></td> <td colspan="7">Nearest to 2.0 + B+F</td> <td></td> <td></td> <td colspan="7">Nearest to 2.0 + B+F</td> </tr> <tr> <td></td> <td colspan="7">Contrast Index</td> <td></td> <td></td> <td colspan="7">Contrast Index</td> </tr> </table>	Film no.	Density Values:								Film no.	Density Values:							3	B + F =		Step no.	Left	Right	Top	Bottom	Average	6	B + F =		Step no.	Left	Right	Top	Bottom	Average		Nearest to 1 + (B+F)									Nearest to 1 + (B+F)								Nearest to 0.5									Nearest to 0.5								Nearest to 2.0 + B+F									Nearest to 2.0 + B+F								Contrast Index									Contrast Index							<table border="1" style="width:100%; border-collapse: collapse;"> <tr> <td style="width:10%;">Film no.</td> <td colspan="7">Density Values:</td> <td style="width:10%;"></td> <td style="width:10%;">Film no.</td> <td colspan="7">Density Values:</td> </tr> <tr> <td>3</td> <td>B + F =</td> <td></td> <td>Step no.</td> <td>Left</td> <td>Right</td> <td>Top</td> <td>Bottom</td> <td>Average</td> <td>6</td> <td>B + F =</td> <td></td> <td>Step no.</td> <td>Left</td> <td>Right</td> <td>Top</td> <td>Bottom</td> <td>Average</td> </tr> <tr> <td></td> <td colspan="7">Nearest to 1 + (B+F)</td> <td></td> <td></td> <td colspan="7">Nearest to 1 + (B+F)</td> </tr> <tr> <td></td> <td colspan="7">Nearest to 0.5</td> <td></td> <td></td> <td colspan="7">Nearest to 0.5</td> </tr> <tr> <td></td> <td colspan="7">Nearest to 2.0 + B+F</td> <td></td> <td></td> <td colspan="7">Nearest to 2.0 + B+F</td> </tr> <tr> <td></td> <td colspan="7">Contrast Index</td> <td></td> <td></td> <td colspan="7">Contrast Index</td> </tr> </table>	Film no.	Density Values:								Film no.	Density Values:							3	B + F =		Step no.	Left	Right	Top	Bottom	Average	6	B + F =		Step no.	Left	Right	Top	Bottom	Average		Nearest to 1 + (B+F)									Nearest to 1 + (B+F)								Nearest to 0.5									Nearest to 0.5								Nearest to 2.0 + B+F									Nearest to 2.0 + B+F								Contrast Index									Contrast Index						
Film no.	Density Values:								Film no.	Density Values:																																																																																																																																																																																																					
3	B + F =		Step no.	Left	Right	Top	Bottom	Average	6	B + F =		Step no.	Left	Right	Top	Bottom	Average																																																																																																																																																																																														
	Nearest to 1 + (B+F)									Nearest to 1 + (B+F)																																																																																																																																																																																																					
	Nearest to 0.5									Nearest to 0.5																																																																																																																																																																																																					
	Nearest to 2.0 + B+F									Nearest to 2.0 + B+F																																																																																																																																																																																																					
	Contrast Index									Contrast Index																																																																																																																																																																																																					
Film no.	Density Values:								Film no.	Density Values:																																																																																																																																																																																																					
3	B + F =		Step no.	Left	Right	Top	Bottom	Average	6	B + F =		Step no.	Left	Right	Top	Bottom	Average																																																																																																																																																																																														
	Nearest to 1 + (B+F)									Nearest to 1 + (B+F)																																																																																																																																																																																																					
	Nearest to 0.5									Nearest to 0.5																																																																																																																																																																																																					
	Nearest to 2.0 + B+F									Nearest to 2.0 + B+F																																																																																																																																																																																																					
	Contrast Index									Contrast Index																																																																																																																																																																																																					

Density Values:									Average of 6 Films									
Film no.	Step no.	1	2	3	4	5	6	Average										
B+ F =																		
Nearest to 1 + (B+F)																		
Nearest to 0.5																		
Nearest to 2.0 + B+F																		
Contrast Index																		

Index Value for Base Fog

Index Value for Speed

Index Value for Contrast

14.3 CROSSOVER WORKSHEET

Processor / Darkroom: Date:.....

Film type:..... Batch / Emulsion No. Date of Expiry:

New Emulsion #					Old Emulsion #				
Film #	Contrast Ind. Low density (LD) Step #	Speed Ind. Mid Density (MD) Step #	Contrast Ind. High density (HD) Step #	Base + Fog (B+F) #	Film #	Contrast Ind. Low density (LD) Step #	Speed Ind. Mid Density (MD) Step #	Contrast Ind. High density (HD) Step #	Base + Fog (B+F) #
1					1				
2					2				
3					3				
4					4				
5					5				
6					6				
Average					1				
Ave. Density Difference: DD = HD - LD =					Ave. Density Difference: DD = HD - LD =				

MD difference between new and old film (New MD - Old MD)	
DD difference between new and old film (New DD - Old DD)	
B + F difference between new and old film (New - Old)	

	MD (Speed Index)	DD (Contrast Index)	B + F
Old operating levels			
Difference between new & old film			
New operating levels			

14.4 CROSSOVER WORKSHEET - EXAMPLE

New Emulsion #					Old Emulsion #				
Film #	Contrast Ind. Low density (LD) Step # 10	Speed Ind. Mid Density (MD) Step # 11	Contrast Ind. High density (HD) Step # 13	Base + Fog (B+F) #	Film #	Contrast Ind. Low density (LD) Step # 10	Speed Ind. Mid Density (MD) Step # 11	Contrast Ind. High density (HD) Step # 13	Base + Fog (B+F) #
1	0.49	1.25	2.39	0.18	1	0.46	1.27	2.33	0.17
2	0.5	1.23	2.43	0.18	2	0.48	1.30	2.30	0.17
3	0.49	1.26	2.40	0.17	3	0.46	1.27	2.28	0.18
4	0.53	1.28	2.41	0.18	4	0.48	1.28	2.32	0.17
5	0.49	1.28	2.43	0.18	5	0.47	1.31	2.31	0.18
6					6				
Average	0.53	1.26	2.41	0.18	1		1.29	2.31	0.17
Ave. Density Difference: DD = HD – LD = 1.91					Ave. Density Difference: DD = HD – LD = 1.84				




MD difference between new and old film (New MD – Old MD)	- 0.03
DD difference between new and old film (New DD – Old DD)	+ 0.07
B + F difference between new and old film (New – Old)	+0.01

	MD (Speed Index)	DD (Contrast Index)	B + F
Old operating levels	1.29	1.84	0.17
Difference between new & old film	- 0.03	+ 0.07	- 0.01
New operating levels	1.26	1.91	0.18

14.5 DAILY PROCESSOR CONTROL CHART:

PROCESSOR: _____ MONTH: _____ YEAR: _____

DENSITOMER MODEL & S/N: _____ FILM TYPE & BATCH NUMBER: _____



DEVELOPER TEMP ° C

2	4	6	8	10	12	14	16	18	20	22	24	26	28	30
---	---	---	---	----	----	----	----	----	----	----	----	----	----	----

FRESH CHEMICALS

2	4	6	8	10	12	14	16	18	20	22	24	26	28	30
---	---	---	---	----	----	----	----	----	----	----	----	----	----	----

14.8 RECORD OF QC TESTS: CASSETTES AND SCREENS

MONTHLY CLEANING AND VISUAL INSPECTION

Cassette size															Comment
Cassette no.															
Date															
Pass/Fail															
Date															
Pass/Fail															
Date															
Pass/Fail															
Date															
Pass/Fail															
Date															
Pass/Fail															
Date															
Pass/Fail															
Date															
Pass/Fail															
Date															
Pass/Fail															
Date															
Pass/Fail															
Date															
Pass/Fail															
Date															

