



Addendum to Media Statement

Please see list of donated equipment below.

No	Description	Qty	Targeted local manufacturers to produce vaccines	Targeted testing	Equipment Manufacturer	Equipment Supplier	Equipment Model
1	High- Performance Liquid Chromatography (HPLC) with variable detectors	1	ASPEN	Potency test	Dionex	Anatech	ICS 6000 Dual system for anionanalysis
2	Optima MAX-XP Tabletop Ultracentrifuge	1	TBI	Potency test	Beckman Coulter	Beckman Coulter	Optima MAX-XP
3	BIO-PLEX 200 SYSTEM WITH HTF	1	ASPEN	Potency test and Identity test	BioRad	Lasec	Bio-Plex 200 system with HTF
4	VITROS® 3600 Immunodiagnostic System	1	ASPEN	Potency test and Identity test	Ortho Clinical Diagnostic	Vacutec	VITROS® 3600 Immunodiagnostic System
5	ELISA plate incubator	2	ASPEN andTBI	Potency test and Identity test	Grant	Monitoring and Control Laboratory	PHMP-4

6	EVOS M7000 Imaging System	1	ASPEN	Potency test and Identity test	Thermo Scientific	Life Technologies	EVOS M7000
7	CO2 Incubator	4	ASPEN	Potency test	Eppendorf	Merck	Cellxpert C170i